PEPTYS Regain
PEPTYS 2
PEPTYS 5

1434 Ed./Rev.: 03 - 20/12/2021

PEPTYS - Instructions for use - Cod. IFU-PEP414

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Eumed S.r.l.
Via Niccolò Tommaseo, 68 - 35131 Padova (ITALIA)

EN

PEPTYS - Solution of low molecular weight peptides (LWPs) derived from hydrolyzed collagen

INSTRUCTIONS FOR USE

PEPTYS - Sterile, single-use device

PEPTYS Regain

PEPTYS 2

PEPTYS 5

DESCRIPTION AND INTENDED USE

Hydrolyzed collagen-derived low molecular weight peptides (LWPs) solution in pre-filled syringe, for infiltrative and injective intra-articular, peri-articular, intra-muscular, intra-dermal and epithelial use.

COMPOSITION

Collagen: low molecular weight peptides (LWPs) derived from hydrolyzed collagen of bovine origin.

Eccipients: vitamin C (stabilized ascorbic acid) and phosphate buffered saline (PBS).

AVAILABLE FORMATS

PEPTYS Regain

- A pre-filled syringe containing 10 ml of LWP solution with a concentration of 1 mg/ml;

PEPTYS 2

- A pre-filled syringe containing 0,5, 1 or 2 ml LWP solution with a concentration of 2 mg/ml;

PEPTYS 5

- A pre-filled syringe containing 0,5, 1 or 2 ml LWP solution with a concentration of 5 mg/ml.

INDICATIONS

PEPTYS is a medical device indicated for the direct reinforcement of the extra-cellular matrix of connective tissues, which is naturally rich in collagen fibers. Its formulation is also recommended for the improvement of articular mobility and for the management of pain symptoms of joints and surrounding connective structures (muscles, tendons, ligaments, menisci) caused by aging, degeneration, trauma, overload, and postural defects. It is also indicated after joint surgery (meniscectomy, ligament reconstruction, articular cartilage repair) to accelerate the functional recovery and to improve post-operative comfort, as well as to counteract the physiological aging and/or the degeneration of connective and epithelial tissues.

MECHANISM OF ACTION

PEPTYS LWPs derive from collagen protein subjected to a hydrolytic process. They can directly support the collagen matrix of connective and epithelial tissues which are compromised because of aging, or damaged due to degenerative and inflammatory pathologies, or after traumatic events.

This mechanical action contributes to the recovery of natural structure and functionality of connective and epithelial tissues, improving for example joint mobility, decreasing symptoms of pain and restoring the

elasticity of extra-cellular matrix.

INSTRUCTIONS FOR USE AND SUGGESTED THERAPEUTIC TREATMENTS

Take the pre-filled syringe out of the sterile envelope. Remove the cap, connect the PEPTYS syringe to a sterile needle with suitable caliber (20-30G, not included in the package) for the specific treatment as suggested here below and proceed to inject under aseptic conditions.

PEPTYS Regain 1 mg/ml

Intra-articular injection following arthroscopy or joint lavage

- Firmly connect the PEPTYS Regain syringe (ref. PEP-110) to a 20-22G needle by turning slightly. Under aseptic conditions, inject into the joint through an arthroscopic portal.

- Suggested treatment: one injection at the end of the arthroscopy surgery.

PEPTYS 2 mg/ml

Intra-articular injection

- Firmly connect the PEPTYS 2 mg/ml syringe (ref. PEP-22 or PEP-21 for small joints) to a 23-28G needle by turning slightly, then proceed with the injection under aseptic conditions.

- Suggested treatment: three injections - day 1, day 15, and day 45.

Intra-tendinous-ligamentous injection

- Firmly connect the PEPTYS 2 mg/ml syringe (ref. PEP-21) to a 26-28G needle by turning slightly, then proceed with the injection under aseptic conditions.

- Suggested treatment: three injections - day 1, day 21, and day 42.

Intra-muscular injection

- Firmly connect the PEPTYS 2 mg/ml syringe (ref. PEP-205 or PEP-21) to a 26-28G needle by turning slightly, then proceed with the injection under aseptic conditions.

- Suggested treatment: two injections - day 1, day 10.

PEPTYS 5 mg/ml

Intra-articular injection

- Firmly connect the PEPTYS 5 mg/ml syringe (ref. PEP-52 or PEP-51 for small joints) to a 23-28G needle by turning slightly, then proceed with the injection under aseptic conditions.

- Suggested treatment: one injection.

Intra-tendinous-ligamentous injection

- Firmly connect the PEPTYS 5 mg/ml syringe (ref. PEP-51) to a 26-28G needle by turning slightly, then proceed with the injection under aseptic conditions.

- Suggested treatment: one injection.

Intra-muscular injection

- Firmly connect the PEPTYS 5 mg/ml syringe (ref. PEP-505 or PEP-51) to a 26-28G needle by turning slightly, then proceed with the injection under aseptic conditions.

- Suggested treatment: one injection.

Intra-dermal and epithelial injection

- Firmly connect the PEPTYS 5 mg/ml syringe (ref. PEP-505, PEP-51 or PEP-52) to a 30G needle by turning slightly, then proceed with the injection under aseptic conditions.

- Suggested treatment: a single injection in one or more sites in the area to be treated.

WARNINGS AND PRECAUTIONS

Read the instructions carefully before the use.

PEPTYS is a Class III Medical Device, for single-use and for individual patient; it cannot be reused, neither resterilized. It must be used exclusively by qualified medical personnel, following proper aseptic procedures.

PEPTYS can be used on patients already treated with injection of hyaluronic acid, platelet-rich plasma (PRP), or other autologous cell concentrates.

It is advisable to take a few hours of rest after the injection and to avoid physical efforts or excessive exercise in the following 24 hours.

In case of intra-articular effusion, it is recommended to drain the joint fluid before injecting PEPTYS.

Do not use PEPTYS if the package is opened or damaged.

Do not inject PEPTYS intravascularly.

Do not use PEPTYS after the expiration date. The expiration date refers to the product unopened and properly stored.

Keep out of the reach of children.

CONTRAINDICATIONS AND SIDE EFFECTS

PEPTYS is biocompatible. It is latex and phthalates free.

PEPTYS must not be used on patients with already known hypersensitivity to bovine collagen or to vitamin C. The product is also contraindicated in the presence of joint or peri-articular infections, hemarthrosis, erythema, or psoriatic spots close to the infiltration area.

Undesired effects referred to LWP injection have not been reported. However, sometimes the injection can provoke a minimal skin reaction or reddening caused by the needle, slight edema, short-term and mild pain.

PRODUCT PACKAGING

Copolyester syringe (CPE), in single or double laminated polyester-aluminium (ALU/OPA-PE) envelope. Information leaflet. The product label is placed on the outer box and on the back of the envelope. The label on the envelope includes 4 peel-off labels for product traceability, to be applied on the patient medical records or other documentation.

STERILIZATION AND STORAGE

PEPTYS is gamma ray sterilized at 25 kGy. The product must be stored in a cool and dry place, at a temperature between 15 and 25 °C (room temperature), avoiding direct exposure to sunlight. If properly stored, the integrity of the package and the sterility of the device are guaranteed for 5 years from the production date (see expiration date on the external label).

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Tiss'You S.r.l.